



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ADCK2 Antibody [#C11266]

**Catalog Number:** C11266

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q7Z695

**Other Names:** Uncharacterized aarF domain-containing protein kinase 2, EC 2.7.11.-, ADCK2, AARF

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ADCK2.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ADCK2 antibody detects endogenous levels of total ADCK2 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

Huang C.Q., Submitted (MAY-2003) to the EMBL/GenBank/DDBJ databases.

Hillier L.W., Nature 424:157-164(2003).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

