



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ATBP3 Antibody [#C11268]

Catalog Number: C11268

Concentration: 1mg/ml

Swiss-Prot No.: Q7Z7A3

Other Names: Cancer-associated gene protein,
ATP-binding domain-containing protein 3, ATPBD3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ATBP3.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ATBP3 antibody detects endogenous levels of total ATBP3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Yousef G.M., Submitted (APR-2003) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Beausoleil S.A., Nat. Biotechnol. 24:1285-1292(2006).

