



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MOL2B Antibody [#C11272]

Catalog Number: C11272

Concentration: 1mg/ml

Swiss-Prot No.: Q86TA1

Other Names: Mps one binder kinase activator-like 2B, Mob1 homolog 2b, Protein Mob3B, MOBKL2B, MOB3B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MOL2B.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MOL2B antibody detects endogenous levels of total MOL2B protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Florindo C.S., Submitted (AUG-2003) to the EMBL/GenBank/DDBJ databases.

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

