



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AKAP14 Antibody [#C11279]

Catalog Number: C11279

Concentration: 1mg/ml

Swiss-Prot No.: Q86UN6

Other Names: A-kinase anchor protein 14, A-kinase anchor protein 14, A-kinase anchor protein 28 kDa, AKAP14, AKAP28

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from interinal of human AKAP14.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AKAP14 antibody detects endogenous levels of total AKAP14 protein.


Reactivity: Human

Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000


References:

Kultgen P.L., Mol. Biol. Cell 13:4156-4166(2002).

The MGC Project Team, Genome Res. 14:2121-2127(2004).



Western blot analysis of extracts from Jurkat cells, K562 cells and HUVEC cells, using AKAP14 antibody.



Immunofluorescence analysis of HeLa cells, using AKAP14 antibody.