



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

GLCTK Antibody [#C11307]

Catalog Number: C11307

Concentration: 1mg/ml

Swiss-Prot No.: Q8IVS8

Other Names: Glycerate kinase, EC 2.7.1.31, HBeAg-binding protein 4, GLYCTK, HBEBP4, LP5910

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GLCTK.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GLCTK antibody detects endogenous levels of total GLCTK protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
IF: 1:100~1:500 ELISA: 1:20000

References:

Guo J.-H., DNA Seq. 17:1-7(2006).


Wan D., Proc. Natl. Acad. Sci. U.S.A. 101:15724-15729(2004).

Ota T., Nat. Genet. 36:40-45(2004).




Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue using GLCTK antibody.



Peptide - +

Western blot analysis of extracts from 3T3 cells, using GLCTK antibody.



Peptide - +

Immunofluorescence analysis of A549 cells, using GLCTK antibody.