



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TUSC5 Antibody [#C11319]

Catalog Number: C11319

Concentration: 1mg/ml

Swiss-Prot No.: Q8IXB3

Other Names: Tumor suppressor candidate 5, Protein located at seventeen-p-thirteen point three 1, TUSC5, LOST1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

TUSC5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: TUSC5 antibody detects endogenous levels of total TUSC5 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:

Konishi H., Oncogene 22:1892-1905(2003).



Western blot analysis of extracts from LOVO cells, using TUSC5 antibody.



Immunofluorescence analysis of HeLa cells, using TUSC5 antibody.