



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TBC1D3 Antibody [#C11330]

Catalog Number: C11330

Concentration: 1mg/ml

Swiss-Prot No.: Q8IZP1

Other Names: TBC1 domain family member 3, Rab GTPase-activating protein PRC17, Prostate cancer gene 17 protein, TRE17 alpha protein, TBC1D3, PRC17, TBC1D3A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

TBC1D3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TBC1D3 antibody detects endogenous levels of total TBC1D3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Pei L., Cancer Res. 62:5420-5424(2002).

Wiemann S., Genome Res. 11:422-435(2001).

Onno M., DNA Cell Biol. 12:107-118(1993).

