



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KIAA1967 Antibody [#C11338]

Catalog Number: C11338

Concentration: 1mg/ml

Swiss-Prot No.: Q8N163

Other Names: Protein KIAA1967, Deleted in breast cancer gene 1 protein, DBC.1, DBC-1, p30 DBC, KIAA1967, DBC1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human KIAA1967.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KIAA1967 antibody detects endogenous levels of total KIAA1967 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Nagase T., DNA Res. 8:319-327(2001).

Hamaguchi M., Proc. Natl. Acad. Sci. U.S.A. 99:13647-13652(2002).

Olsen J.V., Cell 127:635-648(2006).

