



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MAP4K6 Antibody [#C11342]

Catalog Number: C11342

Concentration: 1mg/ml

Swiss-Prot No.: Q8N4C8

Other Names: EC 2.7.11.1, Mitogen-activated protein kinase kinase kinase kinase 6, MAPK/ERK kinase kinase 6, MEK kinase kinase 6, MEKKK 6,

Misshapen/NIK-related kinase, GCK family kinase MiNK

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

MAP4K6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MAP4K6 antibody detects endogenous levels of total MAP4K6 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Dan I., FEBS Lett. 469:19-23(2000).

Beausoleil S.A., Proc. Natl. Acad. Sci. U.S.A. 101:12130-12135(2004).

Olsen J.V., Cell 127:635-648(2006).

