



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## DCLK2 Antibody [#C11345]

**Catalog Number:** C11345

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q8N568

**Other Names:** Serine/threonine-protein kinase DCLK2, EC 2.7.11.1, Doublecortin-like and CAM kinase-like 2, Doublecortin-like kinase 2, DCLK2, DCAMKL2, DCDC3B, DCK2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human DCLK2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** DCLK2 antibody detects endogenous levels of total DCLK2 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Totoki Y., Submitted (MAR-2005) to the EMBL/GenBank/DDBJ databases.

Hillier L.W., Nature 434:724-731(2005).

The German cDNA consortium, Submitted (JUL-2002) to the EMBL/GenBank/DDBJ databases.

