



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SPIC Antibody [#C11346]

Catalog Number: C11346

Concentration: 1mg/ml

Swiss-Prot No.: Q8N5J4

Other Names: Transcription factor Spi-C, SPIC

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human SPIC.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SPIC antibody detects endogenous levels of total SPIC protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Carlsson R., Gene 299:271-278(2002).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

