



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FAKD3 Antibody [#C11055]

Catalog Number: C11055

Concentration: 1mg/ml

Swiss-Prot No.: Q14CZ7

Other Names: FAST kinase domain-containing protein 3, FASTKD3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FAKD3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FAKD3 antibody detects endogenous levels of total FAKD3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Sjoeblom T., Science 314:268-274(2006).

