



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

WTAP Antibody [#C11057]

Catalog Number: C11057

Concentration: 1mg/ml

Swiss-Prot No.: Q15007

Other Names: KIAA0105, Putative pre-mRNA splicing regulator female-lethal (2D) homolog, WT1-associated protein, Wilms' tumor 1-associating protein

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

WTAP.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: WTAP antibody detects endogenous levels of total WTAPv protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Ortega A., J. Biol. Chem. 278:3040-3047(2003).

Nagase T., DNA Res. 2:37-43(1995).

Mungall A.J., Nature 425:805-811(2003).

