



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

BRD3 Antibody [#C11064]

Catalog Number: C11064

Concentration: 1mg/ml

Swiss-Prot No.: Q15059

Other Names: Bromodomain-containing protein 3, KIAA0043, RING3-like protein, RING3L

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human BRD3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: BRD3 antibody detects endogenous levels of total BRD3 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Nomura N., DNA Res. 1:223-229(1994).

Humphray S.J., Nature 429:369-374(2004).

Yu L.-R., J. Proteome Res. 6:4150-4162(2007).

