



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CDK10 Antibody [#C11066]

Catalog Number: C11066

Concentration: 1mg/ml

Swiss-Prot No.: Q15131

Other Names: EC 2.7.11.22, PISSLRE, kinase Cdk10

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CDK10.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CDK10 antibody detects endogenous levels of total CDK10 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Grana X., *Oncogene* 9:2097-2103(1994).

Crawford J., *Genomics* 56:90-97(1999).

Wissing J., *Mol. Cell. Proteomics* 6:537-547(2007).

