



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

IRF4 Antibody [#C11070]

Catalog Number: C11070

Concentration: 1mg/ml

Swiss-Prot No.: Q15306

Other Names: interferon regulatory factor 4; LSIRF; lymphocyte-specific interferon regulatory factor; multiple myeloma oncogene 1; MUM1; NF-EM5v

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human IRF4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: IRF4 antibody detects endogenous levels of total IRF4 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Grossman A., Genomics 37:229-233(1996).

Iida S., Nat. Genet. 17:226-230(1997).

Mungall A.J., Nature 425:805-811(2003).

