



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ANKRD1 Antibody [#C11072]

Catalog Number: C11072

Concentration: 1mg/ml

Swiss-Prot No.: Q15327

Other Names: ALRP, C-193, C193, CARP, CVARP, HA1A2, MCARP, ankyrin repeat domain-containing 1, cardiac ankyrin repeat, cytokine inducible nuclear protein, cytokine-inducible nuclear, liver ankyrin repeat domain 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ANKRD1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ANKRD1 antibody detects endogenous levels of total ANKRD1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Chu W., J. Biol. Chem. 270:10236-10245(1995).

Park J.-H., Cancer Res. 65:2804-2814(2005).

Deloukas P., Nature 429:375-381(2004).

