



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

SURF1 Antibody [#C11080]

Catalog Number: C11080

Concentration: 1mg/ml

Swiss-Prot No.: Q15526

Other Names: Surfeit locus protein 1, SURF1, SURF-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human SURF1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SURF1 antibody detects endogenous levels of total SURF1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Lennard A., DNA Cell Biol. 13:1117-1126(1994).

Zhu Z., Nat. Genet. 20:337-343(1998).

Pequignot M.O., Hum. Mutat. 17:374-381(2001).

