



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SSXT Antibody [#C11081]

Catalog Number: C11081

Concentration: 1mg/ml

Swiss-Prot No.: Q15532

Other Names: SSXT protein, Synovial sarcoma, translocated to X chromosome, SYT protein, SS18, SSXT, SYT

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human

SSXT.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SSXT antibody detects endogenous levels of total SSXT protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Clark J., Nat. Genet. 7:502-508(1994).

Brodin B., Gene 268:173-182(2001).

Fligman I., Am. J. Pathol. 147:1592-1599(1995).

