



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HIPK4 Antibody [#C11365]

Catalog Number: C11365

Concentration: 1mg/ml

Swiss-Prot No.: Q8NE63

Other Names: Homeodomain-interacting protein kinase 4, EC 2.7.11.1, HIPK4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human HIPK4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HIPK4 antibody detects endogenous levels of total HIPK4 protein.


Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
IF: 1:00~1:500 ELISA: 1:20000

References:


Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team, Genome Res. 14:2121-2127(2004).




Peptide - +

Immunohistochemistry analysis of paraffin-embedded human testis tissue using HIPK4 antibody.



Peptide - +

Western blot analysis of extracts from HUVEC cells, using HIPK4 antibody.



Peptide - +

Immunofluorescence analysis of A549 cells, using HIPK4 antibody.