



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ACRBP Antibody [#C11366]

Catalog Number: C11366

Concentration: 1mg/ml

Swiss-Prot No.: Q8NEB7

Other Names: Acrosin-binding protein [Precursor],
Proacrosin-binding protein sp32, Cancer/testis antigen
OY-TES-1, Cancer/testis antigen 23, CT23, ACRBP

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human

ACRBP.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: ACRBP antibody detects endogenous levels
of total ACRBP protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Ono T., Proc. Natl. Acad. Sci. U.S.A. 98:3282-3287(2001).
The MGC Project Team, Genome Res.
14:2121-2127(2004).

