



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RGAG1 Antibody [#C11369]

Catalog Number: C11369

Concentration: 1mg/ml

Swiss-Prot No.: Q8NET4

Other Names: Retrotransposon gag domain-containing protein 1, BJ-HCC-23 tumor antigen, RGAG1, KIAA1318

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human RGAG1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RGAG1 antibody detects endogenous levels of total RGAG1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Xueyuan D., Submitted (JUN-2002) to the EMBL/GenBank/DDBJ databases.

Nagase T., DNA Res. 7:65-73(2000).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

