



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ZNF397 Antibody [#C11371]

Catalog Number: C11371

Concentration: 1mg/ml

Swiss-Prot No.: Q8NF99

Other Names: Zinc finger protein 397, Zinc finger and SCAN domain-containing protein 15, Zinc finger protein 47, ZNF397, ZNF47, ZSCAN15

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human ZNF397.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ZNF397 antibody detects endogenous levels of total ZNF397 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Wu Y., Gene 310:193-201(2003).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Matsuoka S., Science 316:1160-1166(2007).

