



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TAF5 Antibody [#C11082]

Catalog Number: C11082

Concentration: 1mg/ml

Swiss-Prot No.: Q15542

Other Names: Transcription initiation factor TFIID subunit 5, Transcription initiation factor TFIID 100 kDa subunit, TAF(II)100, TAFII-100, TAFII100, TAF5, TAF2D

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TAF5.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TAF5 antibody detects endogenous levels of total TAF5 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:1000

References:

Dubrovskaya V., EMBO J. 15:3702-3712(1996).

Tanese N., Proc. Natl. Acad. Sci. U.S.A. 93:13611-13616(1996).


Tao Y., J. Biol. Chem. 272:6714-6721(1997).



Peptide -

+

Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using TAF5 antibody.



Peptide -

+

Western blot analysis of extracts from A549 cells, using TAF5 antibody.