



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

IPMK Antibody [#C11376]

Catalog Number: C11376

Concentration: 1mg/ml

Swiss-Prot No.: Q8NFU5

Other Names: Inositol polyphosphate multikinase, EC 2.7.1.151, Inositol 1,3,4,6-tetrakisphosphate 5-kinase, IPMK, IMPK

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human IPMK.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: IPMK antibody detects endogenous levels of total IPMK protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Nalaskowski M.M., Biochem. J. 366:549-556(2002).

Chang S.-C., J. Biol. Chem. 277:43836-43843(2002).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

