



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TAF13 Antibody [#C11083]

Catalog Number: C11083

Concentration: 1mg/ml

Swiss-Prot No.: Q15543

Other Names: Transcription initiation factor TFIID subunit 13, Transcription initiation factor TFIID 18 kDa subunit, TAF(II)18, TAFII-18, TAFII18, TAF2K, TAFII18

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TAF13.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TAF13 antibody detects endogenous levels of total TAF13 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:5000

References:

Mengus G., EMBO J. 14:1520-1531(1995).

Gregory S.G., Nature 441:315-321(2006).

Birck C., Cell 94:239-249(1998).



Peptide -

+

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using TAF13 antibody.