



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

ADCK3 Antibody [#C11378]

Catalog Number: C11378

Concentration: 1mg/ml

Swiss-Prot No.: Q8NI60

Other Names: Chaperone activity of bc1 complex-like, mitochondrial [Precursor], Chaperone-ABC1-like, aarF domain-containing protein kinase 3, CABC1, ADCK3, PP265

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ADCK3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADCK3 antibody detects endogenous levels of total ADCK3 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Iiizumi M., Cancer Res. 62:1246-1250(2002).

Wan D., Proc. Natl. Acad. Sci. U.S.A. 101:15724-15729(2004).

Gregory S.G., Nature 441:315-321(2006).

