



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PP14C Antibody [#C11380]

Catalog Number: C11380

Concentration: 1mg/ml

Swiss-Prot No.: Q8TAE6

Other Names: Protein phosphatase 1 regulatory subunit 14C, PKC-potentiated PP1 inhibitory protein, Kinase-enhanced PP1 inhibitor, Serologically defined breast cancer antigen NY-BR-81, PPP1R14C, KEPI

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

PP14C.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PP14C antibody detects endogenous levels of total PP14C protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Scanlan M.J., Cancer Immun. 1:4-4(2001).

Liu Q.-R., J. Biol. Chem. 277:13312-13320(2002).

Mungall A.J., Nature 425:805-811(2003).

