



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

LMO3 Antibody [#C11382]

Catalog Number: C11382

Concentration: 1mg/ml

Swiss-Prot No.: Q8TAP4

Other Names: LIM domain only protein 3,
Neuronal-specific transcription factor DAT1, Rhombotin-3,
LMO3, RBTN3, RBTNL2, RHOM3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human
LMO3.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: LMO3 antibody detects endogenous levels of
total LMO3 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team, Genome Res.
14:2121-2127(2004).

Froni L., J. Mol. Biol. 226:747-761(1992).

