



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR62 Antibody [#G348]

Catalog Number: G348

Concentration: 1mg/ml

Swiss-Prot No.: Q9BZJ7

Other Names: Probable G-protein coupled receptor 62;
hGPCR8; GPR62

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

GPR62.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR62 antibody detects endogenous levels of total GPR62 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Lee D.K., Brain Res. Mol. Brain Res. 86:13-22(2001).

Takeda S., FEBS Lett. 520:97-101(2002).

