



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## GPR82 Antibody [#G354]

**Catalog Number:** G354

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q96P67

**Other Names:** G protein-coupled receptor 8

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human GPR82.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** GPR82 antibody detects endogenous levels of total GPR82 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Lee D.K., Gene 275:83-91(2001).

Takeda S., FEBS Lett. 520:97-101(2002).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

