



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR83 Antibody [#G355]

Catalog Number: G355

Concentration: 1mg/ml

Swiss-Prot No.: Q9NYM4

Other Names: Probable G-protein coupled receptor 83;
G-protein coupled receptor 72; GPR83; GPR72,
KIAA1540

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human GPR83.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR83 antibody detects endogenous levels of total GPR83 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Parker R., Biochim. Biophys. Acta 1491:369-375(2000).

Nagase T., DNA Res. 7:143-150(2000).

Ota T., Nat. Genet. 36:40-45(2004).

