



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

KISSR Antibody [#G373]

Catalog Number: G373

Concentration: 1mg/ml

Swiss-Prot No.: Q969F8

Other Names: KiSS-1 receptor; KiSS-1R; Kisspeptins receptor; Metastin receptor; G-protein coupled receptor 54; Hypogonadotropin-1; hOT7T175; KISS1R; AXOR12, GPR54; OT7T175

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

KISSR.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KISSR antibody detects endogenous levels of total KISSR protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:5000

References:

Ohtaki T., Nature 411:613-617(2001).

Clements M.K., Biochem. Biophys. Res. Commun. 284:1189-1193(2001).

Muir A.I., J. Biol. Chem. 276:28969-28975(2001).

