



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

LGR4 Antibody [#G374]

Catalog Number: G374

Concentration: 1mg/ml

Swiss-Prot No.: Q9BXB1

Other Names: Leucine-rich repeat-containing G-protein coupled receptor 4 [Precursor]; G-protein coupled receptor 48; LGR4; GPR48

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

LGR4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: LGR4 antibody detects endogenous levels of total LGR4 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Loh E.D., Biochem. Biophys. Res. Commun. 282:757-764(2001).

