



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

LPHN1 Antibody [#G377]

Catalog Number: G377

Concentration: 1mg/ml

Swiss-Prot No.: O94910

Other Names: Latrophilin-1; Calcium-independent alpha-latrotoxin receptor 1; Lectomedin-2; LPHN1; KIAA0821, LEC2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human LPHN1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: LPHN1 antibody detects endogenous levels of total LPHN1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:1000

References:

Douangpanya J., Submitted (SEP-2000) to the EMBL/GenBank/DDBJ databases.

Nagase T., DNA Res. 5:355-364(1998).

Brill L.M., Anal. Chem. 76:2763-2772(2004).

