



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

MC4R Antibody [#G384]

Catalog Number: G384

Concentration: 1mg/ml

Swiss-Prot No.: P32245

Other Names: MC4-R; MC4R; melanocortin-4 receptor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MC4R.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MC4R antibody detects endogenous levels of total MC4R protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Gantz I., J. Biol. Chem. 268:15174-15179(1993).

Mountjoy K.G., Mol. Endocrinol. 8:1298-1308(1994).

Hinney A., J. Clin. Endocrinol. Metab. 84:1483-1486(1999).

