



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

MC5R Antibody [#G385]

Catalog Number: G385

Concentration: 1mg/ml

Swiss-Prot No.: P33032

Other Names: MC-2; MC5-R; melanocortin 5 receptor; melanocortin receptor 5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human MC5R.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MC5R antibody detects endogenous levels of total MC5R protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:40000

References:


Chhajlani V., Biochem. Biophys. Res. Commun. 195:866-873(1993).

Griffon N., Biochem. Biophys. Res. Commun. 200:1007-1014(1994).

Fathi Z., Neurochem. Res. 20:107-113(1995).



Western blot analysis of extracts from K562 cells, using MC5R antibody.



Immunofluorescence analysis of MCF-7 cells, using MC5R antibody.