



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

MPRA Antibody [#G388]

Catalog Number: G388

Concentration: 1mg/ml

Swiss-Prot No.: Q86WK9

Other Names: Membrane progesterin receptor alpha; mPR alpha; Progesterin and adipoQ receptor family member VII; PAQR7; MRPA

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MPRA.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MPRA antibody detects endogenous levels of total MPRA protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Zhu Y., Proc. Natl. Acad. Sci. U.S.A. 100:2237-2242(2003).

Tang Y.T., J. Mol. Evol. 61:372-380(2005).

Hammes S.R.; Proc. Natl. Acad. Sci. U.S.A. 100:2168-2170(2003).

