



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## MRGRE Antibody [#G391]

**Catalog Number:** G391

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q86SM8

**Other Names:** Mas-related G-protein coupled receptor member E; G-protein coupled receptor 167; MRGPRES; GPR167, MRGE

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human MRGRE.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** MRGRE antibody detects endogenous levels of total MRGRE protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Vassilatis D.K., Proc. Natl. Acad. Sci. U.S.A. 100:4903-4908(2003).

