



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SMRC2 Antibody [#C11383]

Catalog Number: C11383

Concentration: 1mg/ml

Swiss-Prot No.: Q8TAQ2

Other Names: SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily C member 2, SWI/SNF complex 170 kDa subunit, BRG1-associated factor 170, BAF170, SMARCC2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

SMRC2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SMRC2 antibody detects endogenous levels of total SMRC2 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Wang W., Genes Dev. 10:2117-2130(1996).

Kim J.-E., J. Proteome Res. 4:1339-1346(2005).

Olsen J.V., Cell 127:635-648(2006).

