



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RIN3 Antibody [#C11385]

Catalog Number: C11385

Concentration: 1mg/ml

Swiss-Prot No.: Q8TB24

Other Names: Ras and Rab interactor 3, Ras interaction/interference protein 3, RIN3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RIN3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RIN3 antibody detects endogenous levels of total RIN3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Heilig R., Nature 421:601-607(2003).

Ota T., Nat. Genet. 36:40-45(2004).

Saito K., J. Biol. Chem. 277:3412-3418(2002).

