



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ZNF596 Antibody [#C11388]

Catalog Number: C11388

Concentration: 1mg/ml

Swiss-Prot No.: Q8TC21

Other Names: Zinc finger protein 596, ZNF596

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ZNF596.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: ZNF596 antibody detects endogenous levels of total ZNF596 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100


ELISA: 1:1000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Nusbaum C., Nature 439:331-335(2006).


The MGC Project Team, Genome Res. 14:2121-2127(2004).



Peptide -

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue using ZNF596 antibody.



ZNF596 --

Peptide - +

Western blot analysis of extracts from Jurkat cells, using ZNF596 antibody.