



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AGR3 Antibody [#C11390]

Catalog Number: C11390

Concentration: 1mg/ml

Swiss-Prot No.: Q8TD06

Other Names: AG-3; hAG-3; AG3; Breast cancer membrane protein 11

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human AGR3.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AGR3 antibody detects endogenous levels of total AGR3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:5000

References:

Adam P.J., J. Biol. Chem. 278:6482-6489(2003).


Clark H.F., Genome Res. 13:2265-2270(2003).

Zhang Z., Protein Sci. 13:2819-2824(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using AGR3 antibody.



Peptide - +

Immunofluorescence analysis of COS-7 cells, using AGR3 antibody.