



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ES8L3 Antibody [#C11394]

Catalog Number: C11394

Concentration: 1mg/ml

Swiss-Prot No.: Q8TE67

Other Names: Epidermal growth factor receptor pathway substrate 8-related protein 3; EPS8-like protein 3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ES8L3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ES8L3 antibody detects endogenous levels of total ES8L3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:10000

References:

Tocchetti A., Genomics 81:234-244(2003).


Ota T., Nat. Genet. 36:40-45(2004).

Gregory S.G., Nature 441:315-321(2006).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using ES8L3 antibody.



Peptide - +

Immunofluorescence analysis of HeLa cells, using ES8L3 antibody.