



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

GRK7 Antibody [#C11400]

Catalog Number: C11400

Concentration: 1mg/ml

Swiss-Prot No.: Q8WTQ7

Other Names: G protein-coupled receptor kinase 7 [Precursor], EC 2.7.11.16, G protein-coupled receptor kinase GRK7, GRK7, GPRK7

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

GRK7.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GRK7 antibody detects endogenous levels of total GRK7 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Chen C.-K., Mol. Vis. 7:305-313(2001).

Weiss E.R., J. Neurosci. 21:9175-9184(2001).

