



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ORAV1 Antibody [#C11406]

Catalog Number: C11406

Concentration: 1mg/ml

Swiss-Prot No.: Q8WV07

Other Names: Oral cancer overexpressed protein 1,
Tumor amplified and overexpressed sequence 1,
ORAOV1, TAOS1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ORAV1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ORAV1 antibody detects endogenous levels of total ORAV1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Huang X., Proc. Natl. Acad. Sci. U.S.A. 99:11369-11374(2002).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

