



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CNKR2 Antibody [#C11415]

Catalog Number: C11415

Concentration: 1mg/ml

Swiss-Prot No.: Q8WXI2

Other Names: Connector enhancer of kinase suppressor of ras 2, Connector enhancer of KSR2, CNK2, CNKSR2, CNK2, KIAA0902, KSR2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

CNKR2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CNKR2 antibody detects endogenous levels of total CNKR2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Lanigan T.M., FASEB J. 17:2048-2060(2003).

Nagase T., DNA Res. 5:355-364(1998).

Ross M.T., Nature 434:325-337(2005).

