



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## SEPT1 Antibody [#C11420]

**Catalog Number:** C11420

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q8WYJ6

**Other Names:** Septin-1, LARP, Serologically defined breast cancer antigen NY-BR-24, Peanut-like protein 3, SEPT1, DIFF6, PNUTL3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human SEPT1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** SEPT1 antibody detects endogenous levels of total SEPT1 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Tu Q., Submitted (AUG-1998) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Scanlan M.J., Cancer Immun. 1:4-4(2001).

