



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## CDKL2 Antibody [#C11433]

**Catalog Number:** C11433

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q92772

**Other Names:** Cyclin-dependent kinase-like 2, EC 2.7.11.22, Serine/threonine-protein kinase KKIAMRE, Protein kinase p56 KKIAMRE, CDKL2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

CDKL2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** CDKL2 antibody detects endogenous levels of total CDKL2 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Taglienti C.A., Oncogene 13:2563-2574(1996).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

