



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## CELF-1 Antibody [#C11440]

**Catalog Number:** C11440

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q92879

**Other Names:** CELF-1, Bruno-like protein 2, RNA-binding protein BRUNOL-2, CUG triplet repeat RNA-binding protein 1, CUG-BP1, Deadenylation factor CUG-BP, 50 kDa nuclear polyadenylated RNA-binding protein, Embryo deadenylation element-binding protein homolog, EDEN-BP homolog

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against

synthesized peptide derived from internal of human CELF-1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** CELF-1 antibody detects endogenous levels of total CELF-1 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Timchenko L.T., Nucleic Acids Res. 24:4407-4414(1996).  
Good P.J., J. Biol. Chem. 275:28583-28592(2000).  
Takahashi N., J. Biochem. 130:581-587(2001).

